

Product datasheet for **AR09724PU-L**

PDRG1 (1-133, His-tag) Human Protein


Product data:

Product Type:	Recombinant Proteins
Description:	PDRG1 (1-133, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SSGLVPRGSH</u> MLSPEAERVL RYLVEVEELA EEVLADKRQI VDLDTKRNQN REGLRALQKD LSLSEDVMVC FGNMFIKMPH PETKEMIEKD QDHLDKEIEK LRKQLKVKVN RLFEAQGKPE LKGFNLNPLN QDELKALKVI LKG
Tag:	His-tag
Predicted MW:	17.6 kDa
Concentration:	lot specific
Purity:	>90%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 5 mM DTT, 20% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human PDRG1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_110442</u>
Locus ID:	81572
UniProt ID:	<u>Q9NUG6, A0A384NKN4</u>
Cytogenetics:	20q11.21
Synonyms:	C20orf126; PDRG
Summary:	May play a role in chaperone-mediated protein folding.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:



15% SDS-PAGE (3ug)