

Product datasheet for **AR09733PU-L**

Serum amyloid A-4 protein (SAA4) (21-130, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Serum amyloid A-4 protein (SAA4) (21-130, His-tag) human recombinant protein, 0.25 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH</u> <u>SSGLVPRGSH</u> <u>MWRSFFKEAL</u> QGVGDMGRAY WDIMISNHQN SNRYLYARGN YDAAQRGPGG VWAAKLISRS RVYLQGLIDY YLFGNSSTVL EDSKSNEKAE EWGRSGKDPD RFRPDGLPKK Y
Tag:	His-tag
Predicted MW:	14.9 kDa
Concentration:	lot specific
Purity:	>90%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 5 mM DTT, 0.2M NaCl, 0.5 mM EDTA
Preparation:	Liquid purified protein
Protein Description:	Recombinant human SAA4 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_006503</u>
Locus ID:	6291
UniProt ID:	<u>P35542</u>
Cytogenetics:	11p15.1
Synonyms:	C-SAA; CSAA



[View online »](#)

Summary: Major acute phase reactant. Apolipoprotein of the HDL complex.[UniProtKB/Swiss-Prot Function]

Product images:

