

## Product datasheet for **AR09746PU-N**

### SEPSECS (1-501, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	SEPSECS (1-501, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH</u> <u>SSGLVPRGSH</u> MNRESFAAGE RLVSPAYVRQ GCEARRSHEH LIRLLEKKGK CPENGWDEST LELFLHELAI MDSNNFLGNC GVGEREGRVA SALVARRHYR FIHGIGRSGD ISAVQPKAAG SLLNKITNS LVLDIILAG VHTVANCFVW PMATGMSLTL CFLTLRHKRP KAKYIWPRI DQKSCFKSMI TAGFEPVIE NVLEGDELRT DLKAVEAKVQ ELGPDCILCI HSTTSCFAPR VPDRLEELAV ICANYDIPHI VNNAYGVQSS KCMHLIQQGA RVGRIDAFVQ SLDKNFMVPV GGAIAGFND SFIQEISKMY PGRASASPSL DVLITLLSLG SNGYKLLKE RKEMFSYLSN QIKKLSEAYN ERLHTPHNP ISLAMTLKTL DEHRDKAVTQ LGSMLFTRQV SGARVVPLGS MQTVSGYTFR GFMSHTNNYP CAYLNAASAI GMKMQDVDFL IKRLDRCLKA VRKERSKESD DNYDKTEDVD IEEMALKLDN VLLDTYQDAS S
Tag:	His-tag
Predicted MW:	57.9 kDa
Concentration:	lot specific
Purity:	>90% pure by SDS-PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.5) containing 40% Glycerol, 0.15M NaCl, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant Human SEPSECS protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for 1-2 weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_058651</u>
Locus ID:	51091



[View online »](#)

UniProt ID:	<a href="#">Q9HD40</a>
Cytogenetics:	4p15.2
Synonyms:	LP; PCH2D; SLA; SLA/LP
Summary:	The amino acid selenocysteine is the only amino acid that does not have its own tRNA synthetase. Instead, this amino acid is synthesized on its cognate tRNA in a three step process. The protein encoded by this gene catalyzes the third step in the process, the conversion of O-phosphoseryl-tRNA(Sec) to selenocysteinyI-tRNA(Sec).[provided by RefSeq, Mar 2011]
Protein Pathways:	Aminoacyl-tRNA biosynthesis

**Product images:**