

## Product datasheet for **AR09780PU-N**

### SULT1C2 / SULT1C1 (1-296, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	SULT1C2 / SULT1C1 (1-296, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SSSLVPRGSH</u> MALTSDLGKQ IKLKEVEGTL LQPATVDNWS QIQSFEAKPD DLLICTYPKA GTTWIQEIVD MIEQNGDVEK CQRAIQHRH PFIEWARPPQ PSGVEKAKAM PSPRIKTHL STQLLPPSFV ENNCKFLYVA RNAKDCMVS YHFQRMNHML PDPGTWEEYF ETFINGKVVW GSWFDHVKGW WEMKDRHQIL FLFYEDIKRD PKHEIRKVMQ FMGKKVDETV LDKIVQETSF EKMKENPMTN RSTVSKSILD QSISFMRKG TVGDWKNHFT VAQNERFDEI YRRKMEGTSI NFCMEL
Tag:	His-tag
Predicted MW:	37 kDa
Purity:	>90%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human SULT1C2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_001047</u>
Locus ID:	6819
UniProt ID:	<u>Q00338</u>
Cytogenetics:	2q12.3
Synonyms:	humSULTC2; ST1C1; ST1C2; SULT1C1



[View online »](#)

**Summary:**

Sulfotransferase enzymes catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs, and xenobiotic compounds. These cytosolic enzymes are different in their tissue distributions and substrate specificities. The gene structure (number and length of exons) is similar among family members. This gene encodes a protein that belongs to the SULT1 subfamily, responsible for transferring a sulfo moiety from PAPS to phenol-containing compounds. Two alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]

**Product images:**