

Product datasheet for **AR09830PU-N**

NUDT2 / APAH1 (1-147, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	NUDT2 / APAH1 (1-147, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SSGLVPRGSH</u> MALRACGLII FRRCLIPKVD NNAIEFLLLQ ASDGIHHWTP PKGHVPEGED DLETALRETQ EEAGIEAGQL TIIEGFKREL NYVARNKPKT VIYWLAEVKD YDVEIRLSHE HQAYRWLGLE EACQLAQFKE MKAALQEGHQ FLCSIEA
Tag:	His-tag
Concentration:	lot specific
Purity:	>90%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris buffer (pH 8.0) containing 10% glycerol, 1 mM DTT, 0.1M NaCl
Preparation:	Liquid purified protein
Protein Description:	Recombinant NUDT2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_001152</u>
Locus ID:	318
UniProt ID:	<u>P50583</u>
Cytogenetics:	9p13.3
Synonyms:	APAH1



[View online »](#)

Summary:

This gene encodes a member of the MutT family of nucleotide pyrophosphatases, a subset of the larger NUDIX hydrolase family. The gene product possesses a modification of the MutT sequence motif found in certain nucleotide pyrophosphatases. The enzyme asymmetrically hydrolyzes Ap4A to yield AMP and ATP and is responsible for maintaining the intracellular level of the dinucleotide Ap4A, the function of which has yet to be established. This gene may be a candidate tumor suppressor gene. Alternative splicing has been observed at this locus and four transcript variants, all encoding the same protein, have been identified. [provided by RefSeq, Sep 2011]

Protein Pathways:

Purine metabolism, Pyrimidine metabolism

Product images: