

Product datasheet for **AR09841PU-L**

PPP3R1 / Calcineurin B (1-170, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	PPP3R1 / Calcineurin B (1-170, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Tag:	His-tag
Predicted MW:	21.5 kDa
Concentration:	lot specific
Purity:	>90%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 2mM DTT, 0.1M NaCl
Preparation:	Liquid purified protein
Protein Description:	Recombinant human PPP3R1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_000936
Locus ID:	5534
Cytogenetics:	2p14
Synonyms:	CALNB1; CNB; CNB1
Summary:	Regulatory subunit of calcineurin, a calcium-dependent, calmodulin stimulated protein phosphatase. Confers calcium sensitivity.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Phosphatase



[View online »](#)

Protein Pathways:

Alzheimer's disease, Amyotrophic lateral sclerosis (ALS), Apoptosis, Axon guidance, B cell receptor signaling pathway, Calcium signaling pathway, Long-term potentiation, MAPK signaling pathway, Natural killer cell mediated cytotoxicity, Oocyte meiosis, T cell receptor signaling pathway, VEGF signaling pathway, Wnt signaling pathway

Product images: