

Product datasheet for AR09859PU-S

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com

OriGene Technologies, Inc.

Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

SNRPE (1-92, His-tag) Human Protein

Product data:

Product Type: Recombinant Proteins

Description: SNRPE (1-92, His-tag) human recombinant protein, 10 μg

Species: Human
Expression Host: E. coli

Expression cDNA Clone MGSSHHHHHH SSGLVPRGSH MAYRGQGQKV QKVMVQPINL IFRYLQNRSR IQVWLYEQVN

or AA Sequence: MRIEGCIIGF DEYMNLVLDD AEEIHSKTKS RKQLGRIMLK GDNITLLQSV SN

Tag: His-tag

Predicted MW: 12.9 kDa

Concentration: lot specific

Purity: >95%

Buffer: Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 5 mM DTT, 1 mM EDTA,

30% glycerol

Preparation: Liquid purified protein

Protein Description: Recombinant human SNRPE protein, fused to His-tag at N-terminus, was expressed in E.coli

and purified by using conventional chromatography techniques.

Storage: Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeq: NP 001291393

 Locus ID:
 6635

 UniProt ID:
 P62304

 Cytogenetics:
 1q32.1

Synonyms: HYPT11; Sm-E; SME; snRNP-E





Summary: The protein encoded by this gene is a core component of U small nuclear ribonucleoproteins,

which are key components of the pre-mRNA processing spliceosome. The encoded protein plays a role in the 3' end processing of histone transcripts. This protein is one of the targets in the autoimmune disease systemic lupus erythematosus, and mutations in this gene have been associated with hypotrichosis. Several pseudogenes of this gene have been identified.

[provided by RefSeq, Jun 2016]

Protein Families: Stem cell - Pluripotency

Protein Pathways: Spliceosome

Product images:

