

## Product datasheet for **AR09898PU-L**

### **S100A10 (1-97, His-tag) Human Protein**

#### **Product data:**

Product Type:	Recombinant Proteins
Description:	S100A10 (1-97, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SGLVPRGSH</u> MPSQMEHAME TMMFTFHKFA GDKGYLTKED LRVLMEKEFP GFLENQKDPL AVDKIMKDLQ QCRDGVKGFQ SFFSLIAGLT IACNDYFVVH MKQKGKK
Tag:	His-tag
Predicted MW:	13.3 kDa
Concentration:	lot specific
Purity:	>90%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 1 mM DTT, 40% glycerol, 0.1M NaCl, 0.1 mM PMSF
Preparation:	Liquid purified protein
Protein Description:	Recombinant human S100A10 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_002957</u>
Locus ID:	6281
UniProt ID:	<u>P60903</u>
Cytogenetics:	1q21.3
Synonyms:	42C; ANX2L; ANX2LG; CAL1L; Ca[1]; CLP11; GP11; p10; P11



[View online »](#)

**Summary:**

The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein may function in exocytosis and endocytosis. [provided by RefSeq, Jul 2008]

**Protein Families:**

Transmembrane

**Product images:**