

## Product datasheet for **AR09918PU-N**

### PDCD6 (1-191, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	PDCD6 (1-191, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SSGLVPRGSH</u> MAAYSYPGP GAGPGAAGA ALPDQSFLWN VFQRVDKDRS GVISDTELQQ ALSNGTWTPF NPVTVRSIIS MFDRENKAGV NFSEFTGVWK YITDWQNVFR TYDRDNSGMI DKNELKQALS GFGYRLSDQF HDILIRKFDR QGRGQIAFDD FIQGCIVLQR LTDIFRRYDT DQDGWVQVSY EQYLSMVFSI V
Tag:	His-tag
Concentration:	lot specific
Purity:	>95%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 10 mM sodium citrate (pH 3.5) containing 40% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human PDCD6 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_001254485</u>
Locus ID:	10016
UniProt ID:	<u>O75340, Q53FC3</u>
Cytogenetics:	5p15.33
Synonyms:	ALG-2; ALG2; PEF1B



[View online »](#)

**Summary:**

This gene encodes a calcium-binding protein belonging to the penta-EF-hand protein family. Calcium binding is important for homodimerization and for conformational changes required for binding to other protein partners. This gene product participates in T cell receptor-, Fas-, and glucocorticoid-induced programmed cell death. In mice deficient for this gene product, however, apoptosis was not blocked suggesting this gene product is functionally redundant. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is also located on the short arm of chromosome 5. [provided by RefSeq, May 2012]

**Protein Families:**

Druggable Genome

**Product images:**