

## Product datasheet for **AR09922PU-L**

### **CMPK1 (1-228, His-tag) Human Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	CMPK1 (1-228, His-tag) human recombinant protein, 0.5 mg
<b>Species:</b>	Human
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	<u>MGSSHHHHHH SSGLVPRGSH</u> MLSRCRSGLL HVLGLSFLQ TRRPILLCSP RLMKPLVWFV LGGPGAGKGT QCARIVEKYG YTHLSAGELL RDERKNPDSQ YGELIEKYIK EGKIVPEIT ISLLKREMDQ TMAANAQKNK FLIDGFPRNQ DNLQGWNKTM DGKADVSVL FFDCNNEICI ERCLERKSKS GRSDDNRESL EKRIQTYLQS TKPIIDL YEE MGKVKKIDAS KSVDEVFDEV VQIFDKEG
<b>Tag:</b>	His-tag
<b>Concentration:</b>	lot specific
<b>Purity:</b>	>90%
<b>Buffer:</b>	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 1 mM DTT, 100 mM NaCl
<b>Preparation:</b>	Liquid purified protein
<b>Protein Description:</b>	Recombinant human CMPK1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.
<b>Storage:</b>	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.
<b>RefSeq:</b>	<u>NP_001129612</u>
<b>Locus ID:</b>	51727
<b>UniProt ID:</b>	<u>P30085, A0A494BXC7</u>
<b>Cytogenetics:</b>	1p33
<b>Synonyms:</b>	CK; CMK; CMPK; UMK; UMP-CMPK; UMPK



[View online »](#)

**Summary:** This gene encodes one of the enzymes required for cellular nucleic acid biosynthesis. This enzyme catalyzes the transfer of a phosphate group from ATP to CMP, UMP, or dCMP, to form the corresponding diphosphate nucleotide. Alternate splicing results in both coding and non-coding transcript variants. [provided by RefSeq, Feb 2012]

**Protein Families:** Druggable Genome

**Protein Pathways:** Metabolic pathways, Pyrimidine metabolism

**Product images:**

