

Product datasheet for **AR09926PU-L**

Galectin-7 (1-136, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Galectin-7 (1-136, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SSGLVPRGSH</u> MSNVPHKSSL PEGIRPGTVL RIRGLVPPNA SRFHVNLLCG EEQGSDAALH FNRLDTSEV VFNSKEQGSW GREERGPGVP FQRGQPFEVL IASDDGFKAVVGDAQYHHF RHRLPLARVR LVEVGGDVQL DSVRIF
Tag:	His-tag
Predicted MW:	17.2 kDa
Concentration:	lot specific
Purity:	>95%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human LGALS7 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_002298</u>
Locus ID:	3963
UniProt ID:	<u>P47929</u>
Cytogenetics:	19q13.2
Synonyms:	GAL7; LGALS7A



[View online »](#)

Summary:

The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. Differential and in situ hybridization studies indicate that this lectin is specifically expressed in keratinocytes and found mainly in stratified squamous epithelium. A duplicate copy of this gene (GeneID:653499) is found adjacent to, but on the opposite strand on chromosome 19. [provided by RefSeq, Jul 2008]

Protein Families:

Secreted Protein

Product images: