

Product datasheet for **AR09965PU-L**

Methionine aminopeptidase / MAP (1-264, His-tag) Escherichia coli Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Methionine aminopeptidase / MAP (1-264, His-tag) e. coli recombinant protein, 0.5 mg
Species:	Escherichia coli
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH</u> <u>SSGLVPRGSH</u> MAISIKTPED IEKMRVAGRL AAEVLEMIEP YVKPGVSTGE LDRICNDYIV NEQHAVSACL GYHGYPKSVC ISINEVVCHG IPDDAKLLKD GDIVNIDVTV IKDGFHGDTS KMFIVGKPTI MGERLCRITQ ESLYLALRMV KPGINLREIG AAIQKFVEAE GFSVVREYCG HGIGRGFHEE PQVLHYDSRE TNVVLKPGMT FTIEPMVNAG KKEIRTMKDG WTVKTKDRSL SAQYEHTIVV TDNGCEILTL RKDDTIPAI SHDE
Tag:	His-tag
Predicted MW:	31.5 kDa
Concentration:	lot specific
Purity:	>90% pure by SDS-PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant E.coli map protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Summary:	Methionine aminopeptidases, also designated peptidase M proteins, are members of the M24 family of proteins. This protein removes the amino-terminal methionine residue from nascent polypeptides. The active site of MAP contains two adjacent divalent metal ions connected by a water molecule or hydroxide ion.



[View online »](#)

Product images:

