

Product datasheet for **AR100096L**

Human IgG RapidTiter Strips

Product data:

Product Type: Assay Reagents

Description: Human IgG RapidTiter Strips for rapid detection and semi quantitative measurement of human IgG, and human IgG Fc fused proteins, 100 Strips (4xAR100096)

Stability: Strips should not be used beyond the expiration date. Leave RapidTiter Strips in its canister until just before use. Keep dry. Avoid prolong exposure to humidity.

Product images:




Human IgG RapidTiter Strips



Human IgG RapidTiter Strips
Reference Card

Catalog #: AR100096, AR100096L



Human IgG RapidTiter Strips Reference Card

9620 Medical Center Drive, Rockville, MD 20850 | 888-267-4436 | www.origene.com