

Product datasheet for **AR100097L**

Mouse IgG RapidTiter Strips

Product data:

Product Type:	Assay Reagents
Description:	Mouse IgG RapidTiter Strips for rapid detection and semi quantitative measurement of mouse IgG, and mouse IgG Fc fused proteins, 100 Strips (4xAR100097)
Stability:	Strips should not be used beyond the expiration date. Leave RapidTiter Strips in its canister until just before use. Keep dry. Avoid prolong exposure to humidity.

Product images:






Mouse IgG RapidTiter Strips




Mouse IgG RapidTiter Strips
Reference Card

Catalog #: AR100097, AR100097L



9620 Medical Center Drive, Rockville, MD 20850 | 888-267-4436 | www.origene.com

Mouse IgG RapidTiter Strips Reference Card



9620 Medical Center Drive, Rockville, MD 20850 | 888-267-4436 | www.origene.com

Mouse IgG RapidTiter Strips Reference Card