

Product datasheet for **AR100098**

Rabbit IgG RapidTiter Strips

Product data:

Product Type:	Assay Reagents
Description:	Rabbit IgG RapidTiter Strips for rapid detection and semi quantitative measurement of Rabbit IgG, and Rabbit IgG Fc fused proteins, 25 Strips
Stability:	Strips should not be used beyond the expiration date. Leave RapidTiter Strips in its canister until just before use. Keep dry. Avoid prolong exposure to humidity.

Product images:




Rabbit IgG RapidTiter Strips



**Rabbit IgG RapidTiter Strips
Reference Card**

Catalog #: AR100098, AR100098L



Rabbit IgG RapidTiter Strips Reference Card