

Product datasheet for **AR100099**

Biotin RapidTiter Strips

Product data:

Product Type: Assay Reagents

Description: Biotin RapidTiter Strips for semi-quantitative analysis of biotinylated proteins, 25 Strips

Stability: Strips should not be used beyond the expiration date. Leave RapidTiter Strips in its canister until just before use. Keep dry. Avoid prolong exposure to humidity.

Product images:



Biotin RapidTiter Strips




Biotin RapidTiter Strips Reference Card

Catalog #: AR100099 and AR100099L



Rapid Biotinylation Detection



Biotin RapidTiter Strips Reference Card


9620 Medical Center Drive, Rockville, MD 20850 | 888-267-4436 | www.origene.com

Biotin RapidTiter Strips Reference Card

Catalog #: AR100099 and AR100099L



Semi-quantitative Biotin Labeling Ratio Measurement



Biotin RapidTiter Strips Reference Card

9620 Medical Center Drive, Rockville, MD 20850 | 888-267-4436 | www.origene.com