

Product datasheet for **AR31019PU-N**

Hepatitis B Surface Antigen / HBsAg (ayw) Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Hepatitis B Surface Antigen (HBsAg) (ayw) Recombinant HBsAg subtype ayw. Contains full-length HBsAg (226 amino acids, MW of 25.5kDa)
Expression Host:	Pichia pastoris
Predicted MW:	25.5 kDa
Concentration:	lot specific
Purity:	> 95 % by SDS-PAGE (Coomassie staining). Purified from broken yeast cells using clarification, microfiltration, chromatography, gel filtration, ultra filtration and sterile filtration.
Buffer:	0.2M Sodium Phosphate, pH 7.4 containing 0.03M Sodium Chloride without preservatives
Preparation:	Liquid purified protein
Storage:	Store the antigen at 2-8°C. DO NOT FREEZE! Shelf life: one year from despatch.



[View online »](#)