

Product datasheet for **AR31118PU-N**

DYKDDDDK Epitope Tag Protein

Product data:

Product Type:	Recombinant Proteins
Description:	DYKDDDDK Epitope Tag protein, 5.0 mg
Expression cDNA Clone or AA Sequence:	DYKDDDDK-OH
Predicted MW:	1012.96 Da
Concentration:	N/A
Purity:	>95% by HPLC
Buffer:	State: Lyophilized purified peptide as powder
Preparation:	Purified peptide
Applications:	The DYKDDDDK peptide (Flag) is used as a fusion tag for immunoaffinity chromatography.
Protein Description:	Flag Peptide
Storage:	Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.



[View online »](#)