

Product datasheet for **AR31121PU-N**

Ovalbumin (257-264) -COOH Chicken Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Ovalbumin (257-264) -COOH protein, 5.0 mg
Species:	Chicken
Predicted MW:	963.13 Da
Purity:	>95% by HPLC
Buffer:	State: Purified peptide
Preparation:	Purified peptide
Protein Description:	OVA Peptide (257-264) NH ₂ -SIINFEKL-COOH
Storage:	Store undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_990483
Locus ID:	396058
UniProt ID:	P01012
Synonyms:	SERPINB14
Summary:	ovalbumin (SERPINB14)



[View online »](#)