

## Product datasheet for **AR31122PU-S**

### 6xHistidine (HHHHHH) (-OH) Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	6xHistidine (HHHHHH) -OH protein, 1.0 mg
Expression cDNA Clone or AA Sequence:	HHHHHH-OH (1-letter code) H-His-His-His-His-His-His-OH (3-letter code)
Predicted MW:	840.91 Da
Purity:	>95% by HPLC
Buffer:	State: Purified peptide
Preparation:	Purified peptide
Protein Description:	HexHis
Storage:	Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.



[View online »](#)