

# Product datasheet for AR39019PU-N

#### 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

## **UBE2W (1-151, His-tag) Human Protein**

**Product data:** 

**Product Type: Recombinant Proteins** 

**Description:** UBE2W (1-151, His-tag) human recombinant protein, 0.1 mg

Species: Human E. coli **Expression Host:** 

**Expression cDNA Clone** 

MGSSHHHHHH SSGLVPRGSH MASMQKRLQK ELLALQNDPP PGMTLNEKSV QNSITQWIVD MEGAPGTLYE GEKFQLLFKF SSRYPFDSPQ VMFTGENIPV HPHVYSNGHI CLSILTEDWS or AA Sequence:

PALSVQSVCL SIISMLSSCK EKRRPPDNSF YVRTCNKNPK KTKWWYHDDT C

Tag: His-tag Predicted MW: 19.5 kDa Concentration: lot specific **Purity:** >90%

**Buffer:** Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 2 mM DTT, 0.1M NaCl

**Preparation:** Liquid purified protein

**Protein Description:** Recombinant human UBE2W protein, fused to His-tag at N-terminus, was expressed in E.coli

and purified by using conventional chromatography techniques.

Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Storage:

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeq: NP 001001481

Locus ID: 55284

**UniProt ID:** Q96B02, X6REH9

Cytogenetics: 8q21.11

Synonyms: UBC-16; UBC16



#### UBE2W (1-151, His-tag) Human Protein - AR39019PU-N

Summary: This gene encodes a nuclear-localized ubiquitin-conjugating enzyme (E2) that, along with

ubiquitin-activating (E1) and ligating (E3) enzymes, coordinates the addition of a ubiquitin moiety to existing proteins. The encoded protein promotes the ubiquitination of Fanconi anemia complementation group proteins and may be important in the repair of DNA damage. There is a pseudogene for this gene on chromosome 1. Alternative splicing results in

multiple transcript variants. [provided by RefSeq, Aug 2012]

**Protein Families:** Transcription Factors

**Protein Pathways:** Ubiquitin mediated proteolysis

## **Product images:**

