

Product datasheet for **AR39027PU-L**

DYNLL2 (1-89, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	DYNLL2 (1-89, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SSGLVPRGSH</u> MSDRKAVIKN ADMSEDMQQD AVDCATQAME KYNIEKDIAA YIKKEFDKKY NPTWHCIVGR NFGSYVTHET KHFIYFYLQ VAILLFKSG
Tag:	His-tag
Predicted MW:	12.5 kDa
Concentration:	lot specific
Purity:	>90%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 1 mM DTT, 30% glycerol, 0.2M NaCl
Preparation:	Liquid purified protein
Protein Description:	Recombinant human DYNLL2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_542408</u>
Locus ID:	140735
UniProt ID:	<u>Q96FJ2</u>
Cytogenetics:	17q22
Synonyms:	Dlc2; DNCL1B; RSPH22



[View online »](#)

Summary:

Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be involved in linking dynein to cargos and to adapter proteins that regulate dynein function. Cytoplasmic dynein 1 acts as a motor for the intracellular retrograde motility of vesicles and organelles along microtubules. May play a role in changing or maintaining the spatial distribution of cytoskeletal structures (By similarity).
[UniProtKB/Swiss-Prot Function]

Product images: