

Product datasheet for **AR39036PU-N**

AGR3 (22-166, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	AGR3 (22-166, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SSGLVPRGSH MGSMIAIKKE KRPPQTL SRG WGDDITWVQT YEEGLFYAQK SKKPLMVIHH LEDCQYSQAL KKVFAQNEEI QEMAQNKFIM LNLMHETTDK NLSPDGQYVP RIMFVDPSLT VRADIAGRYS NRLTYEPRD LPLLIENMKK ALRLIQSEL</u>
Tag:	His-tag
Predicted MW:	19.5 kDa
Concentration:	lot specific
Purity:	>90%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human AGR3 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_789783</u>
Locus ID:	155465
UniProt ID:	<u>Q8TD06</u>
Cytogenetics:	7p21.1
Synonyms:	AG-3; AG3; BCMP11; hAG-3; HAG3; PDIA18



[View online »](#)

Summary:

This gene encodes a member of the disulfide isomerase (PDI) family of endoplasmic reticulum (ER) proteins that catalyze protein folding and thiol-disulfide interchange reactions. The encoded protein has an N-terminal ER-signal sequence, a catalytically active thioredoxin domain, and a C-terminal ER-retention sequence. This gene is expressed in ciliated airway epithelial cells and, in mouse, plays a role in ciliary beat frequency in multiciliated cells. This gene is also over-expressed in breast, ovarian, and prostate cancers. [provided by RefSeq, Dec 2016]

Protein Families:

Secreted Protein

Product images: