

Product datasheet for **AR39070PU-L**

ARL11 / ARLTS1 (1-196, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	ARL11 / ARLTS1 (1-196, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH SSGLVPRGSH</u> MGSVNSRGHK AEAQVMMGL DSAGKTTLLY KLKGHQLVET LPTVGFNVEP LKAPGHVSLT LWDVGGQAPL RASWKDYLEG TDILVYVLDS TDEARLPESA AELTEVLNDP NMAGVPFLVL ANKQEAPDAL PLLKIRNRLS LERFQDHCWE LRGCSALTGE GLPEALQSLW SLLKSRSCMC LQARAHGAER GDSKRS
Tag:	His-tag
Predicted MW:	23.6 kDa
Concentration:	lot specific
Purity:	>90%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 5 mM DTT, 30% glycerol, 100 mM NaCl, 1 mM EDTA
Preparation:	Liquid purified protein
Protein Description:	Recombinant human ARL11 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_612459</u>
Locus ID:	115761
UniProt ID:	<u>Q969Q4</u>
Cytogenetics:	13q14.2
Synonyms:	ARLTS1



[View online »](#)

Summary:

This gene encodes a tumor suppressor related to the ADP-ribosylation factor (ARF) family of proteins. The encoded protein may play a role in apoptosis in a caspase-dependent manner. Polymorphisms in this gene have been associated with some familial cancers. [provided by RefSeq, May 2010]

Product images: