

Product datasheet for **AR39079PU-L**

Glucagon (90-180, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Glucagon (90-180, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH</u> <u>SSGLVPRGSH</u> <u>MKRHDEFERH</u> AEGTFTSDVS SYLEGQAAKE FIAWLVKGRG RRDFPEEVAI VEELGRRHAD GSFSDDEMNTI LDNLAARDFI NWLIQTKITD RK
Tag:	His-tag
Predicted MW:	12.8 kDa
Concentration:	lot specific
Purity:	>95%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 1 mM DTT, 50% glycerol, 0.2M NaCl
Preparation:	Liquid purified protein
Protein Description:	Recombinant human Glucagon protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_002045</u>
Locus ID:	2641
UniProt ID:	<u>P01275</u>
Cytogenetics:	2q24.2
Synonyms:	GLP-1; GLP1; GLP2; GRPP



[View online »](#)

Summary:

The protein encoded by this gene is actually a preproprotein that is cleaved into four distinct mature peptides. One of these, glucagon, is a pancreatic hormone that counteracts the glucose-lowering action of insulin by stimulating glycogenolysis and gluconeogenesis. Glucagon is a ligand for a specific G-protein linked receptor whose signalling pathway controls cell proliferation. Two of the other peptides are secreted from gut endocrine cells and promote nutrient absorption through distinct mechanisms. Finally, the fourth peptide is similar to glicentin, an active enteroglucagon. [provided by RefSeq, Jul 2008]

Protein Families:

Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein

Product images: