

Product datasheet for **AR39111PU-L**

PCNP (1-178, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	PCNP (1-178, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH</u> <u>SSGLVPRGSH</u> <u>MGSHMADGKA</u> GDEKPEKSQR AGAAGGPEEE AEKPVKTKTV SSSNGGESSS RSAEKRSAEE EAADLPTKPT KISKFGFAIG SQTTKKASAI SIKLGGSSKPK ETVPTLAPKT LSVAAAFNED EDSEPEEMPP EAKMRMKNIG RDTPTSAGPN SFNKGKHGFS DNQKLWERNI KSHLGNVHDQ DN
Tag:	His-tag
Predicted MW:	21.5 kDa
Concentration:	lot specific
Purity:	>85%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.5) containing 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human PCNP protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<u>NP_001307324</u>
Locus ID:	57092
UniProt ID:	<u>Q8WW12</u>
Cytogenetics:	3q12.3
Summary:	May be involved in cell cycle regulation.[UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

