

## Product datasheet for **AR39113PU-N**

### Ribonuclease HI / rnhA (1-155, His-tag) Escherichia coli Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Ribonuclease HI / rnhA (1-155, His-tag) e. coli recombinant protein, 0.1 mg
Species:	Escherichia coli
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	<u>MGSSHHHHHH</u> SSGLVPRGSH MGSMLKQVEI FTDGSCLGNP GPGGYGAILR YRGREKTFSA GYTRTTNNRM ELMAAIVALE ALKEHCEVIL STDSQYVRQG ITQWIHNWKK RGWKTADKKP VKNVDLWQRL DAALGQHQIK WEWVKGHAGH PENERCDELA RAAAMNPTLE DTGYQVEV
Tag:	His-tag
Concentration:	lot specific
Purity:	>95%
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 2 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant e.coli rnhA protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Summary:	rnhA is an endonuclease that specifically degrades the RNA of RNA-DNA hybrids. Localized to the nucleus, this protein mediates the removal of Okazaki fragment RNA primers that are present on the lagging strand during DNA replication. rnhA catalyzes the endonucleolytic cleavage of RNA to a 5'-phosphomonoester and is able to bind magnesium or manganese as cofactors.



[View online »](#)

Product images:

