

Product datasheet for AR39138PU-N

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com

OriGene Technologies, Inc.

EU: info-de@origene.com CN: techsupport@origene.cn

DERA (1-318, His-tag) Human Protein

Product data:

Product Type: Recombinant Proteins

Description: DERA (1-318, His-tag) human recombinant protein, 50 μg

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

MGSSHHHHHH SSGLVPRGSH MSAHNRGTEL DLSWISKIQV NHPAVLRRAE QIQARRTVKK EWQAAWLLKA VTFIDLTTLS GDDTSSNIQR LCYKAKYPIR EDLLKALNMH DKGITTAAVC VYPARVCDAV KALKAAGCNI PVASVAAGFP AGQTHLKTRL EEIRLAVEDG ATEIDVVINR SLVLTGQWEA LYDEIRQFRK ACGEAHLKTI LATGELGTLT NVYKASMIAM MAGSDFIKTS

SLVLTGQWEA LYDEIRQFRK ACGEAHLKTI LATGELGTLT NVYKASMIAM MAGSDFIKTS TGKETVNATF PVAIVMLRAI RDFFWKTGNK IGFKPAGGIR SAKDSLAWLS LVKEELGDEW

LKPELFRIGA STLLSDIERQ IYHHVTGRYA AYHDLPMS

Tag:His-tagPredicted MW:37.3 kDaConcentration:lot specific

Purity: >90%

Buffer: Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl, 1 mM DTT

Preparation: Liquid purified protein

Protein Description: Recombinant human DERA protein, fused to His-tag at N-terminus, was expressed in E.coli

and purified by using conventional chromatography.

Storage: Store undiluted at 2-8°C for up to two weeks or (in aliquots) at -20°C or -70°C for longer.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeq: NP 001287708

Locus ID:51071UniProt ID:E9PPM8Cytogenetics:12p12.3

Synonyms: CGI-26; DEOC





Summary: Catalyzes a reversible aldol reaction between acetaldehyde and D-glyceraldehyde 3-

phosphate to generate 2-deoxy-D-ribose 5-phosphate. Participates in stress granule (SG) assembly. May allow ATP production from extracellular deoxyinosine in conditions of energy

deprivation.[UniProtKB/Swiss-Prot Function]

Protein Families: Druggable Genome

Protein Pathways: Pentose phosphate pathway

Product images:

