

Product datasheet for **AR50014PU-N**

HSP10 (1-102) (recombinant) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	HSP10 (1-102) (recombinant) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MAGQAFRKFL PLFDRVLVER SAAETVTKGG IMLPEKSQ GK VLQATWAVG SGSKGKGGEI QPVSVKVGD K VLLPEYGGTK VLLDDKDYFL FRDGDILGKY VD
Predicted MW:	10,9 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris pH 8.0, 50 mM NaCl
Endotoxin:	< 1.0 EU per 1 microgram of protein (determined by LAL method)
Preparation:	Liquid purified protein
Protein Description:	Recombinant human HSP10 was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_002148
Locus ID:	3336
UniProt ID:	P61604 , A0A384N6A4
Cytogenetics:	2q33.1
Synonyms:	CPN10; EPF; GROES; HSP10



[View online »](#)

Summary:

This gene encodes a major heat shock protein which functions as a chaperonin. Its structure consists of a heptameric ring which binds to another heat shock protein in order to form a symmetric, functional heterodimer which enhances protein folding in an ATP-dependent manner. This gene and its co-chaperonin, HSPD1, are arranged in a head-to-head orientation on chromosome 2. Naturally occurring read-through transcription occurs between this locus and the neighboring locus MOBKL3.[provided by RefSeq, Feb 2011]

Protein Families:

Druggable Genome, Stem cell - Pluripotency

Product images: