

Product datasheet for **AR50047PU-S**

ARL3 (1-182, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	ARL3 (1-182, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSHMGLLSI LRKLKSAPDQ EVRILLGLD NAGKTTLLKQ LASEDISHIT PTQGFNIKSV QSQGFKLNWV DIGGQRKIRP YWKNYFENTD ILIYVIDSAD RKRFEETGQE LAELLEEEKL SCVPVLIFAN KQDLLTAAPA SEIAEGLNLH TIRDRVWQIQ SCSALTGEGV QDGMNWCKN VNAKKK
Tag:	His-tag
Predicted MW:	23.0 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 30% glycerol, 0.1M NaCl, 1 mM DTT, 0.1 mM PMSF
Preparation:	Liquid purified protein
Protein Description:	Recombinant human ARL3 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_004302
Locus ID:	403
UniProt ID:	P36405
Cytogenetics:	10q24.32
Synonyms:	ARFL3; JBTS35; RP83



[View online »](#)

Summary: ADP-ribosylation factor-like 3 is a member of the ADP-ribosylation factor family of GTP-binding proteins. ARL3 binds guanine nucleotides but lacks ADP-ribosylation factor activity. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Product images:

