

Product datasheet for **AR50061PU-S**

Angiopoietin-1 (20-498) Human Protein

Product data:

| | |
|--|---|
| Product Type: | Recombinant Proteins |
| Description: | Angiopoietin-1 (20-498) human recombinant protein, 0.1 mg |
| Species: | Human |
| Expression Host: | E. coli |
| Expression cDNA Clone or AA Sequence: | MSNQRRSPEN SGRRYNRIQH GQCAYTFILP EHDGNCREST TDQYNTNALQ RDPAPHVEPDF SSQKLQHLEH VMENYTQWLQ KLENYIVENM KSEMAQIQN AVQNHTATML EIGTSLLSQT AEQTRKLTDV ETQVLNQTSR LEIQLENSL STYKLEKQLL QQTNEILKIHK EKNLSLEHKI LEMEGKHKEE LDTLKEEKEN LQGLVTRQTY IIQELEKQLN RATTNNSVLQ KQQLLEMDTV HNLVNLCTKE GVLLKGGKRE EEKPFRCAD VYQAGFNKSG IYTIYINNMP EPKKVFCNMD VNGGGWTVIQ HREDGSLDFQ RGWKEYKMGF GNPSGEYWLG NEFIFAITSQ RQYMLRIELM DWEGNRAYSQ YDRFHIGNEK QNYRLYLKGH TGTAGKQSSL ILHGADFSTK DADNDNCMCK CALMLTGGWW FDACGPSNLN GMFYTAGQNH GKLNGIKWHY FKGPSYSLRS TTMMIRPLDF |
| Predicted MW: | 55.6 kDa |
| Concentration: | lot specific |
| Purity: | >85% by SDS - PAGE |
| Buffer: | Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) 1M Urea, 5% Glycerol. |
| Preparation: | Liquid purified protein |
| Protein Description: | Recombinant human ANGPT1 protein was expressed in E.coli and purified by using conventional chromatography techniques. |
| Storage: | Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: one year from despatch. |
| RefSeq: | NP_001137 |
| Locus ID: | 284 |
| UniProt ID: | Q15389 |
| Cytogenetics: | 8q23.1 |



[View online »](#)

Synonyms: AGP1; AGPT; ANG1; HAE5

Summary: This gene encodes a secreted glycoprotein that belongs to the angiotensin family. Members of this family play important roles in vascular development and angiogenesis. All angiotensins bind with similar affinity to an endothelial cell-specific tyrosine-protein kinase receptor. The protein encoded by this gene is a secreted glycoprotein that activates the receptor by inducing its tyrosine phosphorylation. It plays a critical role in mediating reciprocal interactions between the endothelium and surrounding matrix and mesenchyme and inhibits endothelial permeability. The protein also contributes to blood vessel maturation and stability, and may be involved in early development of the heart. Mutations in this gene are associated with hereditary angioedema. [provided by RefSeq, Aug 2020]

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein

Product images:

