

Product datasheet for **AR50095PU-S**

UBE2D3 (1-149, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	UBE2D3 (1-149, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MLSNRKCLSK ELSDLARDPP AQCSAGPVGDMFMHWQATIM GPNDSFYQGG VFFLTIHFPT DYPFKPPKVA FTTRIYHPNI NSNGSICLDI LRSQWSPALT ISKVLLSICS LLCDPNPDDP LVPEIARIYK TDRDKYNRIS REWTQKYAM
Tag:	His-tag
Predicted MW:	19.1 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 40% glycerol, 0.15M NaCl, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human UBE2D3 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001287724
Locus ID:	7323
UniProt ID:	D6RAH7
Cytogenetics:	4q24
Synonyms:	E2(17)KB3; UBC4/5; UBCH5C



[View online »](#)

Summary:

The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme functions in the ubiquitination of the tumor-suppressor protein p53, which is induced by an E3 ubiquitin-protein ligase. [provided by RefSeq, Jan 2017]

Protein Pathways:

Ubiquitin mediated proteolysis

Product images: