

## Product datasheet for **AR50131PU-N**

### PSMB4 / PROS26 (46-264, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	PSMB4 / PROS26 (46-264, His-tag) human recombinant protein, 0.25 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MTQNPMVTGT SVLGVKFEGG VWIADMLGS YGSLARFRNI SRIMRVNNT MLGASGDYAD FQYLKQVLGQ MVIDEELLGD GHSYSPRAIH SWLTRAMYSR RSKMNPWNT MVIGGYADGE SFLGYVDMLG VAYEAPSLAT GYGAYLAQPL LREVLEKQPV LSQTEARDLV ERCMRVLYYR DARSYNRFQI ATVTEKGVEI EGPLSTETNW DIAHMISGFE
Tag:	His-tag
Predicted MW:	26.6 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 1 mM DTT, 30% glycerol, 0.1M NaCl
Preparation:	Liquid purified protein
Protein Description:	Recombinant human PSMB4 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_002787</a>
Locus ID:	5692
UniProt ID:	<a href="#">P28070</a> , <a href="#">A0A140VK46</a>
Cytogenetics:	1q21.3
Synonyms:	HN3; HsN3; PRAAS3; PROS-26; PROS26



[View online »](#)

**Summary:**

The proteasome is a multicatalytic proteinase complex with a highly ordered ring-shaped 20S core structure. The core structure is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. An essential function of a modified proteasome, the immunoproteasome, is the processing of class I MHC peptides. This gene encodes a member of the proteasome B-type family, also known as the T1B family, that is a 20S core beta subunit. [provided by RefSeq, Jul 2008]

**Protein Families:**

Druggable Genome, Protease

**Protein Pathways:**

Proteasome

**Product images:**