

Product datasheet for **AR50169PU-S**

G3BP2 (1-449, His-tag) Human Protein

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | G3BP2 (1-449, His-tag) human recombinant protein, 0.1 mg |
| Species: | Human |
| Expression Host: | E. coli |
| Expression cDNA Clone or AA Sequence: | MGSSHHHHHH SSGLVPRGSH MGSHEMMEKP SPLLVGREFV RQYYTLLNKA PEYLHRYFGR NSSYVHGGVD ASGKQPQEAIV GQNDIHHKVL SLNFSECHTK IRHVDAAHATL SDGVVVQVMG LLNSGQPER KFMQTFVLAP EGSVPNKFYV HNDMFRYEDE VFGDSEPELD EESEDEVEEE QEERQSPSEP VQENANSQYV EAHPTVNGIE EPLEESSHEP EPEPESETKT EELKPQVEEK NLEEELEEKST TPPPAEPVSL PQEPPKPRVE AKPEVQSQPP RVREQRPRER PGFPPRGPRP GRGDMEQNDN DNRRIIRYPD SHQLFVGNLP HDIDENELKE FFMSFGNVVE LRINTKGVGG KLPNFGFVVF DDSEPVQRIL IAKPIMFRGE VRLNVEEKKT RAARERETRG GGDDRRDIRR NDRGPGGPRG IVGGGMMRDR DGRGPPPRGG MAQKLGSGRG TGQMEGRFTG QRR |
| Tag: | His-tag |
| Predicted MW: | 53.3 kDa |
| Concentration: | lot specific |
| Purity: | >90% by SDS - PAGE |
| Buffer: | Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 30% glycerol, 0.2M NaCl, 2 mM DTT, 0.1 mM PMSF |
| Preparation: | Liquid purified protein |
| Protein Description: | Recombinant human G3BP2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography. |
| Storage: | Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: one year from despatch. |
| RefSeq: | NP_036429 |
| Locus ID: | 9908 |



[View online »](#)

UniProt ID: [Q9UN86](#), [A0A024RDE5](#)

Cytogenetics: 4q21.1

Summary: Scaffold protein that plays an essential role in cytoplasmic stress granule formation which acts as a platform for antiviral signaling.[UniProtKB/Swiss-Prot Function]

Protein Families: Druggable Genome

Product images:

