

Product datasheet for **AR50178PU-S**

PTPMT1 (28-201, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	PTPMT1 (28-201, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSHEMKVPR AHRDWYHRID PTVLLGALPL RSLTRQLVQD ENVRGVITMN EYETRFLCN SSQEWKRLGV EQLRLSTVDM TGIPTLDNLQ KGVQFALKYQ SLGQCVYVHC KAGRSRSATM VAAyliQVHK WSPEEAVRAI AKIRSYIHIR PGQLDVLKEF HKQITARATK DGTfVIsKT
Tag:	His-tag
Predicted MW:	22.6 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1 mM DTT, 0.15M NaCl
Preparation:	Liquid purified protein
Protein Description:	Recombinant human PTPMT1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001137456
Locus ID:	114971
UniProt ID:	Q8WUK0
Cytogenetics:	11p11.2
Synonyms:	DUSP23; MOSP; PLIP; PNAS-129



[View online »](#)

Summary:

Lipid phosphatase which dephosphorylates phosphatidylglycerophosphate (PGP) to phosphatidylglycerol (PG) (By similarity). PGP is an essential intermediate in the biosynthetic pathway of cardiolipin, a mitochondrial-specific phospholipid regulating the membrane integrity and activities of the organelle (By similarity). Has also been shown to display phosphatase activity toward phosphoprotein substrates, specifically mediates dephosphorylation of mitochondrial proteins, thereby playing an essential role in ATP production (By similarity). Has probably a preference for proteins phosphorylated on Ser and/or Thr residues compared to proteins phosphorylated on Tyr residues (By similarity). Probably involved in regulation of insulin secretion in pancreatic beta cells (By similarity). May prevent intrinsic apoptosis, probably by regulating mitochondrial membrane integrity (PubMed:24709986).[UniProtKB/Swiss-Prot Function]

Protein Families:

Druggable Genome

Product images: