

Product datasheet for **AR50267PU-S**

RPL8 (1-257, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	RPL8 (1-257, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGRVIRGQRK GAGSVFRAHV KHRKGAARLR AVDFAERHGY IKGIVKDIH DPGRGAPLAK VFRDPYRFK KRTELFIAAE GIHTGQFVYC GKKAQLNIGN VLPVGTMP TIVCCLEEKP GDRGKLARAS GNYATVISHN PETKKTRVKL PSGSKKVISS ANRAVGVVA GGGRIDKPIL KAGRAYHKYK AKRNCWPRVR GVAMNPVEHP FGGGNHQHIG KPSTIRRDAP AGRKVLIAA RRTGRLRGTK TVQEKEN
Tag:	His-tag
Predicted MW:	30.2 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 2 mM DTT, 50% glycerol, 300 mM NaCl, 2 mM EDTA
Preparation:	Liquid purified protein
Protein Description:	Recombinant human RPL8 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by conventional chromatography, after refolding of the isolated inclusion bodies in a renaturation buffer.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_000964
Locus ID:	6132
UniProt ID:	P62917
Cytogenetics:	8q24.3



[View online »](#)

Synonyms: L8

Summary: Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L2P family of ribosomal proteins. It is located in the cytoplasm. In rat, the protein associates with the 5.8S rRNA, very likely participates in the binding of aminoacyl-tRNA, and is a constituent of the elongation factor 2-binding site at the ribosomal subunit interface. Alternatively spliced transcript variants encoding the same protein exist. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2008]

Protein Pathways: Ribosome

Product images:

