

Product datasheet for **AR50343PU-N**

NKIRAS2 (1-191, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	NKIRAS2 (1-191, His-tag) human recombinant protein, 0.25 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSHEMGKSCK VVCGQASVG KTSILEQLLY GNHVVGSEMI ETQEDIYVGS IETDRGVREQ VRFYDTRGLR DGAELPRHCF SCTDGYVLVY STDSRESFQR VELLKKEIDK SKDKKEVTIV VLGKCDLQE QRRVDPDVAQ HWAKSEKVKL WEVSVADRRS LLEPFVYLAS KMTQPQSKSA FPLSRKNKGS GSLDG
Tag:	His-tag
Predicted MW:	24.0 kDa
Concentration:	lot specific
Purity:	90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 1 mM DTT, 10% glycerol, 0.1M NaCl
Preparation:	Liquid purified protein
Protein Description:	Recombinant human NKIRAS2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001001349
Locus ID:	28511
UniProt ID:	Q9NYR9 , A0A024R1Z4
Cytogenetics:	17q21.2
Synonyms:	kappaB-Ras2; KBRAS2



[View online »](#)

Summary:

Atypical Ras-like protein that acts as a potent regulator of NF-kappa-B activity by preventing the degradation of NF-kappa-B inhibitor beta (NFKBIB) by most signals, explaining why NFKBIB is more resistant to degradation. May act by blocking phosphorylation of NFKBIB and nuclear localization of p65/RELA NF-kappa-B subunit. It is unclear whether it acts as a GTPase. Both GTP- and GDP-bound forms block phosphorylation of NFKBIB (By similarity). [UniProtKB/Swiss-Prot Function]

Product images: