

## Product datasheet for **AR50382PU-S**

### FABP12 (1-140, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	FABP12 (1-140, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGS HMIDQLQ GTWKSISCEN SEDYMKELGI GRASRKLGR LAKPTVTISTD GDVITIKTKS IFKNNEISFK LGEEFEEITP GGHKTKSKVT LDKESLIQVQ DWDGKETTIT RKLVDGKMV ESTVNSVICT RTYEKVSNS VSNS
Tag:	His-tag
Predicted MW:	18 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 7.5) containing 0.15M NaCl, 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human FABP12 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_001098751</a>
Locus ID:	646486
UniProt ID:	<a href="#">A6NFH5</a>
Cytogenetics:	8q21.13
Summary:	May play a role in lipid transport.[UniProtKB/Swiss-Prot Function]



[View online »](#)

## Product images:

