

Product datasheet for AR50418PU-S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Vacuolar-sorting protein SNF8 (1-258, His-tag) Human Protein

Product data:

Product Type: Recombinant Proteins

Description: Vacuolar-sorting protein SNF8 (1-258, His-tag) human recombinant protein, 0.1 mg

Species: Human E. coli **Expression Host:**

Expression cDNA Clone

MGSSHHHHHH SSGLVPRGSH MGSHMHRRGV GAGAIAKKKL AEAKYKERGT VLAEDQLAQM or AA Sequence:

SKQLDMFKTN LEEFASKHKQ EIRKNPEFRV QFQDMCATIG VDPLASGKGF WSEMLGVGDF

YYELGVQIIE VCLALKHRNG GLITLEELHQ QVLKGRGKFA QDVSQDDLIR AIKKLKALGT GFGIIPVGGT

YLIQSVPAEL NMDHTVVLQL AEKNGYVTVS EIKASLKWET ERARQVLEHL LKEGLAWLDL

QAPGEAHYWL PALFTDLYSQ EITAEEAREA LP

Tag: His-tag Predicted MW: 31.4 kDa Concentration: lot specific

Purity: >90% by SDS - PAGE

Buffer: Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer, pH 8.0, 50% glycerol, 2 mM DTT, 200 mM NaCl

Liquid purified protein Preparation:

Protein Description: Recombinant human SNF8 protein, fused to His-tag at N-terminus, was expressed in E.coli

and purified by using conventional chromatography techniques.

Storage: Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeq: NP 001304121

Locus ID: 11267 **UniProt ID:** Q96H20

Cytogenetics: 17q21.32

Synonyms: Dot3; EAP30; VPS22



Vacuolar-sorting protein SNF8 (1-258, His-tag) Human Protein – AR50418PU-S

Summary: The protein encoded by this gene is a component of the endosomal sorting complex required

for transport II (ESCRT-II), which regulates the movement of ubiquitinylated transmembrane proteins to the lysosome for degradation. This complex also interacts with the RNA

polymerase II elongation factor (ELL) to overcome the repressive effects of ELL on RNA polymerase II activity. Several transcript variants encoding different isoforms have been

found for this gene. [provided by RefSeq, Nov 2015]

Protein Families: Transcription Factors

Protein Pathways: Endocytosis

Product images:

