

Product datasheet for **AR50422PU-N**

ANXA8L2 (1-327, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	ANXA8L2 (1-327, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSHEMAWWKA WIEQEGVTVK SSSHFNPPDP AETLYKAMKG IGTNEQAIID VLTKRSNTQR QQIAKSFKAQ FGKDLTETLK SELSGKFERL IVALMYPPYR YEAKELHDAM KGLGTKEGVI IEILASRTKN QLREIMKAYE EDYSSLEED IQADTSGYLE RILVCLLQGS RDDVSSFVDP ALALQDAQDL YAAGENIRGT DEMKFITILC TRSATHLLRV FEEYEKIANK SIEDSIKSET HGSLEEAMLT VVKCTQNLHS YFAERLYYAM KGAGTRDGTL IRNIVSRSEI DLNLIKCHFK KMYGKTLSSM IMEDTSGDYK NALLSLVGSD P
Tag:	His-tag
Predicted MW:	39.4 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human ANXA8L2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Locus ID:	728113
Cytogenetics:	10q11.22
Synonyms:	ANXA8; ANXA8L2; bA145E20.2; VAC-beta



[View online »](#)

Summary:

This gene encodes a member of the annexin family of evolutionarily conserved Ca²⁺ and phospholipid binding proteins. The encoded protein may function as an anticoagulant that indirectly inhibits the thromboplastin-specific complex. Overexpression of this gene has been associated with acute myelocytic leukemia. A highly similar duplicated copy of this gene is found in close proximity on the long arm of chromosome 10. [provided by RefSeq, Apr 2014]

Product images: