

Product datasheet for **AR50438PU-S**

APIP (1-242, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	APIP (1-242, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMSGCDA REGDCCSRRC GAQDKEHPRY LIPELCKQFY HLGWWTGTGG GISLKHGDEI YIAPSGVQKE RIQPEDMFVC DINEKDISGP SPSKLLKKSQ CTPLFMNAYT MRGAGAVIHT HSKAAVMATL LFPGREFKIT HQEMIKGIKK CTSGGYRYD DMLVPIIEN TPEEKDLKDR MAHAMNEYPD SCAVLVRRHG VYVWGETWEK AKTMCECYDY LFDIAVSMKK VGLDPSQLPV GENIV
Tag:	His-tag
Predicted MW:	29.7 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human APIP protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_057041
Locus ID:	51074
UniProt ID:	Q96GX9
Cytogenetics:	11p13
Synonyms:	APIP2; CGI-29; CGI29; hAPIP; MMRP19



[View online »](#)

Summary: APIP is an APAF1 (MIM 602233)-interacting protein that acts as a negative regulator of ischemic/hypoxic injury (Cho et al., 2004 [PubMed 15262985]).[supplied by OMIM, Dec 2008]

Protein Pathways: Cysteine and methionine metabolism, Metabolic pathways

Product images:

