

Product datasheet for **AR50471PU-S**

Tropomodulin-3 (TMOD3) (1-352, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Tropomodulin-3 (TMOD3) (1-352, His-tag) human recombinant protein, 50 µg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSHPALPFR KDLEKYKDL D EDELLGNLSE TELKQLETVL DDLDPENALL PAGFRQKNQT SKSTTGPFDR EHLLSYLEKE ALEHKDREY VPYTGEKKGK IFIPKQKPVQ TFTEEKVSLD PELEEALSA SDELCDLAA ILGMHNLITN TKFCNIMGSS NGVDQEHFSN VVKGEKILPV FDEPPNPTNV EESLKRKEN DAHLVEVNLN NIKNIPIPTL KDFAKALETN THVKCFSLAA TRSNDPVATA FAEMLKVNKT LKSLNVESNF ITGVGILALI DALRDNETLA ELKIDNQRQQ LGTAVELEMA KMLEENTNIL KFGYQFTQQG PRTRAANAIT KNNDLVRKRR VEGDHQ
Tag:	His-tag
Predicted MW:	42 kDa
Concentration:	lot specific
Purity:	>80% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human TMOD3 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_055362
Locus ID:	29766
UniProt ID:	Q9NYL9
Cytogenetics:	15q21.2
Synonyms:	Ubiquitous tropomodulin, U-Tmod



[View online »](#)

Summary:

Blocks the elongation and depolymerization of the actin filaments at the pointed end. The Tmod/TM complex contributes to the formation of the short actin protofilament, which in turn defines the geometry of the membrane skeleton (By similarity).[UniProtKB/Swiss-Prot Function]

Product images: