

## Product datasheet for **AR50519PU-S**

### Arfaptin-2 (1-341, His-tag) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Arfaptin-2 (1-341, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMTDGILG KAATMEIPIH GNGEARQLPE DDGLEQDLQQ VMVSGPNLNE TSIVSGGYGG SGDGLIPTGS GRHPSHSTTP SGPGDEVARG IAGEKFDIVK KWGINTYKCT KQLLSERFGR GSRTVDLELE LQIELLRETK RKYESVLQLG RALTAHLYSL LQTQHALGDA FADLSQKSPE LQEFGYNAE TQKLLCKNGE TLLGAVNFFV SSINTLVTKT MEDTLMTVKQ YEAARLEYDA YRTDLEELSL GPRDAGTRGR LESAQATFQA HRDKYEKLRG DVAIKLKLE ENKIKVMHKQ LLLFHNAVSA YFAGNQKQLE QTLQQFNIKL RPPGAEKPSW LEEQ
Tag:	His-tag
Predicted MW:	40.2 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human ARFIP2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	<a href="#">NP_001229783</a>
Locus ID:	23647
UniProt ID:	<a href="#">A0A087X1E4</a> , <a href="#">B4DUZ3</a>
Cytogenetics:	11p15.4
Synonyms:	POR1



[View online »](#)

**Summary:** Putative target protein of ADP-ribosylation factor. Involved in membrane ruffling.  
[UniProtKB/Swiss-Prot Function]

**Protein Families:** Druggable Genome

**Product images:**

