

Product datasheet for AR50543PU-S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ACTA2 / aortic smooth muscle Actin (3-377, His-tag) Human Protein

Product data:

Product Type: Recombinant Proteins

Description: ACTA2 / aortic smooth muscle Actin (3-377, His-tag) human recombinant protein, 50 μg

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

 $\frac{\text{MGSSHHHHHH}}{\text{SSGLVPRGSH}} \frac{\text{MGSHM}}{\text{EEEDS}} \text{TALVCDNGSG} \text{ LCKAGFAGDD APRAVFPSIV} \\ \text{GRPRHQGVMV GMGQKDSYVG DEAQSKRGIL TLKYPIEHGI ITNWDDMEKI WHHSFYNELR} \\$

VAPEEHPTLL TEAPLNPKAN REKMTQIMFE TFNVPAMYVA IQAVLSLYAS GRTTGIVLDS

GDGVTHNVPI YEGYALPHAI MRLDLAGRDL TDYLMKILTE RGYSFVTTAE REIVRDIKEK LCYVALDFEN EMATAASSSS LEKSYELPDG QVITIGNERF RCPETLFQPS FIGMESAGIH ETTYNSIMKC DIDIRKDLYA

NNVLSGGTTM YPGIADRMQK EITALAPSTM KIKIIAPPER KYSVWIGGSI LASLSTFQQM

WISKQEYDEA GPSIVHRKCF

Tag: His-tag
Predicted MW: 44.4 kDa
Concentration: lot specific

Purity: >85% by SDS - PAGE

Buffer: Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol

Preparation: Liquid purified protein

Protein Description: Recombinant Human ACTA2 protein, fused to His-tag at N-terminus, was expressed in *E.coli*

and purified by using conventional chromatography techniques.

Storage: Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeg: NP 001135417

Locus ID: 59

UniProt ID: <u>P62736</u>, <u>D2JYH4</u>

Cytogenetics: 10q23.31





Synonyms: ACTSA

Summary: This gene encodes one of six different actin proteins. Actins are highly conserved proteins

that are involved in cell motility, structure, integrity, and intercellular signaling. The encoded protein is a smooth muscle actin that is involved in vascular contractility and blood pressure homeostasis. Mutations in this gene cause a variety of vascular diseases, such as thoracic aortic disease, coronary artery disease, stroke, and Moyamoya disease, as well as

aortic disease, coronary artery disease, stroke, and Moyamoya disease, as well as multisystemic smooth muscle dysfunction syndrome. [provided by RefSeq, Sep 2017]

Protein Pathways: Vascular smooth muscle contraction

Product images:

