

Product datasheet for **AR50619PU-N**

TRAPPC2 (1-140, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	TRAPPC2 (1-140, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSMSGSFYF VIVGHHDNPV FEMEFLPAGK AESKDDHRHL NQFIAHAALD LVDENMWLSN NMYLKTVDKF NEWFVSAFVT AGHMRFIMLH DIRQEDGIKN FFTDVYDLYI KFSMNPYEP NSPIRSSAFD RKVQFLGKKH LLS
Tag:	His-tag
Predicted MW:	18.8 kDa
Concentration:	lot specific
Purity:	>90% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human TRAPPC2 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001011658
Locus ID:	6399
UniProt ID:	P0DI81 , P0DI82 , Q6IBE5
Cytogenetics:	Xp22.2
Synonyms:	hYP38334; MIP2A; SEDL; SEDT; TRAPPC2P1; TRS20; ZNF547L



[View online »](#)

Summary:

The protein encoded by this gene is thought to be part of a large multi-subunit complex involved in the targeting and fusion of endoplasmic reticulum-to-Golgi transport vesicles with their acceptor compartment. In addition, the encoded protein can bind c-myc promoter-binding protein 1 and block its transcriptional repression capability. Mutations in this gene are a cause of spondyloepiphyseal dysplasia tarda (SED). A processed pseudogene of this gene is located on chromosome 19, and other pseudogenes are found on chromosomes 8 and Y. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Mar 2010]

Protein Families:

Druggable Genome, Transcription Factors

Product images: