

## Product datasheet for AR50627PU-S

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## AP3 complex subunit sigma-1 / AP3S1 (1-193, His-tag) Human Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** AP3 complex subunit sigma-1 / AP3S1 (1-193, His-tag) human recombinant protein, 0.1 mg

Species: Human
Expression Host: E. coli

Expression cDNA Clone MGSSHHHHHH SSGLVPRGSH MIKAILIFNN HGKPRLSKFY QPYSEDTQQQ IIRETFHLVS

or AA Sequence: KRDENVCNFL EGGLLIGGSD NKLIYRHYAT LYFVFCVDSS ESELGILDLI QVFVETLDKC FENVCELDLI

FHVDKVHNIL AEMVMGGMVL ETNMNEIVTQ IDAQNKLEKS EAGLAGAPAR AVSAVKNMNL

PEIPRNINIG DISIKVPNLP SFK

Tag: His-tag
Predicted MW: 23.8 kDa
Concentration: lot specific

Purity: >90% by SDS - PAGE

**Buffer:** Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol, 1 mM

DTT

**Preparation:** Liquid purified protein

**Protein Description:** Recombinant human AP3S1 protein, fused to His-tag at N-terminus, was expressed in E.coli

and purified by using conventional chromatography techniques.

**Storage:** Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer.

Avoid repeated freezing and thawing.

**Stability:** Shelf life: one year from despatch.

**RefSeq:** NP 001002924

**Locus ID:** 1176

UniProt ID: Q92572

Cytogenetics: 5q22.3-q23.1

Synonyms: CLAPS3; Sigma3A





Summary: This gene encodes a subunit of the AP3 adaptor complex. This complex functions in the

formation of subcellular vesicles budded from the Golgi body. Several related pseudogenes of


this gene have been found. Alternative splicing results in multiple transcript variants

encoding different isoforms. [provided by RefSeq, Dec 2015]

**Protein Families:** Druggable Genome

**Protein Pathways:** Lysosome

## **Product images:**



15% SDS-PAGE (3ug)