

Product datasheet for AR50644PU-N

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com

OriGene Technologies, Inc.

EU: info-de@origene.com CN: techsupport@origene.cn

RAET1E (31-225, His-tag) Human Protein

Product data:

Product Type: Recombinant Proteins

Description: RAET1E (31-225, His-tag) human recombinant protein, 0.5 mg

Species: Human E. coli **Expression Host:**

Expression cDNA Clone

MGSSHHHHHH SSGLVPRGSH MGSMHSLCFN FTIKSLSRPG QPWCEAQVFL NKNLFLQYNS or AA Sequence: DNNMVKPLGL LGKKVYATST WGELTQTLGE VGRDLRMLLC DIKPQIKTSD PSTLQVEMFC

QREAERCTGA SWQFATNGEK SLLFDAMNMT WTVINHEASK IKETWKKDRG LEKYFRKLSK

GDCDHWLREF LGHWEAMPEP TVSPVNASDI HWSSSSLPD

Tag: His-tag Predicted MW: 24.9 kDa **Concentration:** lot specific

>90% by SDS - PAGE **Purity:**

Buffer: Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.4M UREA, 10% glycerol

Preparation: Liquid purified protein

Protein Description: Recombinant human RAET1E protein, fused to His-tag at N-terminus, was expressed in E.coli.

Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Storage:

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeq: NP 001230254

Locus ID: 135250 **UniProt ID:** Q8TD07 **Cytogenetics:** 6q25.1

Synonyms: bA350J20.7; LETAL; N2DL-4; NKG2DL4; RAET1E2; RL-4; ULBP4






Summary:

This gene belong to the RAET1 family, which consists of major histocompatibility complex (MHC) class I-related genes located in a cluster on chromosome 6q24.2-q25.3. This and RAET1G protein differ from other RAET1 proteins in that they have type I membrane-spanning sequences at their C termini rather than glycosylphosphatidylinositol anchor sequences. This protein functions as a ligand for NKG2D receptor, which is expressed on the surface of several types of immune cells, and is involved in innate and adaptive immune responses. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Aug 2011]

Protein Families: Druggable Genome, Transmembrane
Protein Pathways: Natural killer cell mediated cytotoxicity

Product images:



15% SDS-PAGE (3ug)-