

## Product datasheet for AR50677PU-N

## 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com

OriGene Technologies, Inc.

EU: info-de@origene.com CN: techsupport@origene.cn

## CD3 (23-126, His-tag) Human Protein

**Product data:** 

**Product Type: Recombinant Proteins** 

Description: CD3 (23-126, His-tag) human recombinant protein, 0.25 mg

Species: Human E. coli **Expression Host:** 

**Expression cDNA Clone** 

MGSSHHHHHH SSGLVPRGSH MGSDGNEEMG GITQTPYKVS ISGTTVILTC PQYPGSEILW QHNDKNIGGD EDDKNIGSDE DHLSLKEFSE LEQSGYYVCY PRGSKPEDAN FYLYLRARVC or AA Sequence:

**ENCMEMD** 

His-tag Tag: Predicted MW: 14.1 kDa Concentration: lot specific

**Purity:** >90% by SDS - PAGE

**Buffer:** Presentation State: Purified

State: Liquid purified protein

Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 10% glycerol

Preparation: Liquid purified protein

**Protein Description:** Recombinant human CD3E protein, fused to His-tag at N-terminus, was expressed in E.coli

and purified by using conventional chromatography techniques.

Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Storage:

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

RefSeq: NP 000724

Locus ID: 916 **UniProt ID:** P07766 Cytogenetics: 11q23.3

Synonyms: IMD18; T3E; TCRE





**Summary:** 

The protein encoded by this gene is the CD3-epsilon polypeptide, which together with CD3-gamma, -delta and -zeta, and the T-cell receptor alpha/beta and gamma/delta heterodimers, forms the T-cell receptor-CD3 complex. This complex plays an important role in coupling antigen recognition to several intracellular signal-transduction pathways. The genes encoding the epsilon, gamma and delta polypeptides are located in the same cluster on chromosome 11. The epsilon polypeptide plays an essential role in T-cell development. Defects in this gene cause immunodeficiency. This gene has also been linked to a susceptibility to type I diabetes in women. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome, Transmembrane

**Protein Pathways:** Hematopoietic cell lineage, Primary immunodeficiency, T cell receptor signaling pathway

## **Product images:**

