

Product datasheet for **AR50719PU-N**

AICDA / AID (1-198, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	AICDA / AID (1-198, His-tag) human recombinant protein, 0.5 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MDSLLMNRK FLYQFKNVRW AKGRRETYLC YVVKRRDSAT SFSLDFGYLR NKNNGCHVELL FLRYISDWDL DPGRCYRVTW FTSWSPCYDC ARHVADFLRG NPNLRLRIFT ARLYFCEDRK AEPEGLRRLH RAGVQIAIMT FKDYFYCWNT FVENHERTFK AWEGLHENS V RLSRQLRRIL LPLYEVDDL R DAFRTLGL
Tag:	His-tag
Predicted MW:	26.1 kDa
Concentration:	lot specific
Purity:	>85% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol
Preparation:	Liquid purified protein
Protein Description:	Recombinant human AICDA protein, fused to His-tag at N-terminus, was expressed in E.coli.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001317272
Locus ID:	57379
UniProt ID:	Q9GZX7 , Q7Z599
Cytogenetics:	12p13.31
Synonyms:	AID; ARP2; CDA2; HEL-S-284; HIGM2



[View online »](#)

Summary:


This gene encodes a RNA-editing deaminase that is a member of the cytidine deaminase family. AICDA is specifically expressed and active in germinal center-like B cells. In the germinal center, AICDA is involved in somatic hypermutation, gene conversion, and class-switch recombination of immunoglobulin genes. An epigenetic role in neoplastic transformation and lymphoma progression has been experimentally ascribed to AICDA using mouse models. Defects in this gene are the cause of autosomal recessive hyper-IgM immunodeficiency syndrome type 2 (HIGM2). [provided by RefSeq, Jul 2020]

Protein Families:

Druggable Genome

Protein Pathways:

Primary immunodeficiency

Product images:

15% SDS-PAGE (3ug)