

Product datasheet for **AR50735PU-S**

BOLA1 (21-137, His-tag) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	BOLA1 (21-137, His-tag) human recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	MGSSHHHHHH SSGLVPRGSH MGSQGSAGSG AIGPVEAAIR TKLEEALSPE VLELRNESGG HAVPPGSETH FRVAVVSSRF EGLSPLQRHR LVHAALAEEL GGPVHALAIQ ARTPAQWREN SQLDTSPCL GGNKKTGT
Tag:	His-tag
Predicted MW:	14 kDa
Concentration:	lot specific
Purity:	>95% by SDS - PAGE
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 20 mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 20% glycerol, 1 mM DTT
Preparation:	Liquid purified protein
Protein Description:	Recombinant human BOLA1 protein, fused to His-tag at N-terminus, was expressed in E.coli and purified by using conventional chromatography techniques.
Storage:	Store undiluted at 2-8°C for one week or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001307954
Locus ID:	51027
UniProt ID:	Q9Y3E2
Cytogenetics:	1q21.2
Synonyms:	CGI-143



[View online »](#)

Summary:

Acts as a mitochondrial iron-sulfur (Fe-S) cluster assembly factor that facilitates (Fe-S) cluster insertion into a subset of mitochondrial proteins (By similarity). Probably acts together with the monothiol glutaredoxin GLRX5 (PubMed:27532772). May protect cells against oxidative stress (PubMed:22746225).[UniProtKB/Swiss-Prot Function]

Protein Families:

Transcription Factors

Product images: